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
for compressed gas driven tools

Contents

Use.....	6
Design.....	6
Function and operation.....	6
Technical data.....	7
Important information.....	7
Maintenance / Service.....	7
Assembly information.....	8
Installation dimensions rail system.....	10
Installation dimensions flush-mounting / hollow wall..	10
Installation dimensions ceiling support.....	10
Spare parts.....	11

Airmotor-Outlets for compressed gas driven tools

Use

to drive compressed air driven implements in enclosed areas; spent air is passed into the open air via an exhaust air line

Design

Take-offs for compressed air driven implements are connected to a separate piping network for compressed air (700 to 1000 kPa). Connection of the take-off to the implement is achieved via a coaxial hose. The inner hose passes compressed air to the implement, the outer hose passes return air to the exhaust air duct

Function and operation

Only connectors for the Air Motor can be coupled into the take-off.

Push the connector into the take-off until an audible click is heard. Coupling in the connector opens a valve. Compressed air flows through the supply line to the implement. The compressed air is used here and spent air fed back to the plug-in coupling. A line takes it from there to the exhaust air duct.

Airmotor-Outlets

for compressed gas driven tools

Technical data

Type:	All-metal constructions, stainless steel housing with front panel, one-handed coupling and uncoupling, release via female connector
Connections:	pressure gas line: Cu-tube 8x1 exhaust air line: Cu-tube 15x1
Capacity:	at least 350 l/min at compressed air pressure of 700 to 1000 kPa
Color:	black / white designation
Installation:	- flush-mounting or hollow wall mounting - ceiling supply systems - wall supply systems dimensions see annex
Guidelines:	DIN EN 737-3 DIN 13260-3 (draft)

Important information

- Take-offs and connectors to DIN 13260 T3!
- Never use tools when connecting the unit and the accessories; all connections can be tightened or loosened by hand.
- Use only original GREGGERSEN spare parts

Service/Maintenance

External disinfection with a cloth is sufficient for all chrome-plated brass parts. Clean plastic parts with a common disinfectant solution and follow manufacturers instructions. Let outlets and plug dry completely before using them again.

Prior to each use, make a visual check of the outlets (external damages). Outlets and plug should be controlled at least once a year (capacity, function); in very frequently used areas accordingly more frequently!


Failures and damages may only be remedied by authorized expert personnel. For repairs, use only original GREGGERSEN spare parts!

Airmotor-Outlets

for compressed gas driven tools

Assembly notes

- The flange has two fixing bores (pitch circle 42.5 mm) for installing the take-off (2) in a housing.
- The bezel (5) is fixed with two screws (6). The gas type plate (7) is glued to the bezel.
- The front panel (3) is fixed to the housing with two screws (16).
- The 6 mm compressed air hose (14) is attached (exception: fix compressed air connection, ceiling supply unit, with clamp). The 13 mm exhaust air hose (15) must be fixed with two clamps.

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GREGGERSEN Gasetechnik GmbH